**REF** 161001, 161006, 161009,  
161010, 1610010, 1610060,  
1610090, 1610100

English | 10107351/01  
Software version 3.0.x

**Store this Addendum together with your ventilator Operator's Manual.**

This Addendum summarizes changes and corrections to the *HAMILTON-C1/T1/MR1 Operator's Manuals*.

## 1 Corrections to the HAMILTON-C1/T1/MR1 Operator's Manuals

### 1.1 Using the COM1 RS-232 port

*HAMILTON-C1/T1 only.* RS-232 connection over the COM1 port is only available when using the communication Y-cable (PN 10077038).

### 1.2 Alarm troubleshooting

*HAMILTON-C1/T1 only.* The following alarm is *not* described in the ventilator *Operator's Manual*.

Table 1. Alarms and other messages

Alarm	Definition
CO2: Poor signal	<i>Low priority.</i> The CO2 sensor signal quality is poor.

**Action needed**

- Check patient condition.
- Check CO2 sensor and adapter connections.
- Ensure that airway adapters are not in a horizontal position relative to the floor to reduce accumulation of patient secretions. If accumulation occurs, remove the adapter, rinse with water, and reconnect.

### 1.3 About closed-suctioning maneuvers

When performing a closed-suctioning maneuver, ventilation continues and the current settings do *not* need to be adjusted.

### 1.4 Setting up a bluetooth connection

*HAMILTON-C1/T1 only.* When setting up the connection of a mobile device using the bluetooth connection type, the text "Pairing available" is *not* displayed in the System > Settings > Connectivity > Bluetooth window.

### 1.5 Vt/Weight in Controls window

Due to space constraints in some languages, the Vt/Weight parameter is shown as Vt/kg in the lower left corner of the Controls window, when either the (S)CMV+ /APVcmv or (S)IMV+/APVsimv mode is selected (Neonatal patient group only).

## 1.6 HAMILTON-C1/T1 water protection

*HAMILTON-C1/T1 only.* The water protection rating of the HAMILTON-C1 and HAMILTON-T1 ventilators has changed, and the device labels have been updated accordingly.

The following table describes the HAMILTON-C1 and HAMILTON-T1 ingress protection ratings.

### Symbols used on device labels and packaging

Device	Symbol	Definition
HAMILTON-C1	<b>IP22</b>	Protected against dripping water when the device is tilted to a maximum of 15 degrees, and from solid particles larger than 12.5 mm.
HAMILTON-T1	<b>IP54</b>	Protected from water spray from any direction, and from limited dust ingress.

## 1.7 HAMILTON-T1 environmental requirements

*HAMILTON-T1 only.* When the device is *not* stored in its original packaging, the permitted storage temperature range is -15°C to 60°C.

## 1.8 HAMILTON-T1 standards and approvals

*HAMILTON-T1 only.* Table 15-20 of the ventilator Operator's Manual has been updated as follows:

Table 2. Standards

IEC 60601-1	Medical electrical equipment, Part 1: General requirements for basic safety and essential performance. The device classification is: Class II, Type BF applied part (ventilator breathing system, VBS, CO2 sensor including CO2 module connector, and SpO2 sensor including SpO2 adapter), continuous operation.
IEC 60601-1-2	Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance <ul style="list-style-type: none"> <li>• Collateral standard: Electromagnetic disturbances</li> <li>• Requirements and tests</li> </ul>
IEC 60601-1-10	Medical electrical equipment - Part 1-10: General requirements for basic safety and essential performance Collateral standard: Requirements for the development of physiologic closed-loop controllers
IEC 60601-1-12	Medical electrical equipment - Part 1-12: General requirements for basic safety and essential performance Collateral standard: Requirements for medical electrical equipment and medical electrical systems intended for use in the emergency medical services environment

ISO 80601-2-12	Medical electrical equipment - Part 2-12: Particular requirements for the basic safety and essential performance of critical care ventilators
CAN/CSA-C22.2 No. 60601.1	Medical electrical equipment: General requirements for safety
ANSI/AAMI ES 60601-1	Medical electrical equipment - Part 1: General requirements for basic safety and essential performance
EN ISO 60601-1	Anaesthetic and respiratory equipment - conical connectors - Part 1: Cones and sockets
EN ISO 5359	Low-pressure hose assemblies for use with medical gases
EN ISO 80601-2-55	Medical electrical equipment - Part 2-55: Particular requirements for the basic safety and essential performance of respiratory gas monitors
MIL-STD-461F	Control of electromagnetic interference
MIL-STD-810G	Low pressure (altitude)
EN 1789	Medical vehicles and their equipment - Road ambulances
EN 794-3	Lung ventilators - Part 3: Particular requirements for emergency and transport ventilators

## 2 Corrections to the Pulse Oximetry Instructions for Use for HAMILTON-C1/T1

### 2.1 Operating conditions

The stated operating conditions of the SpO2 sensor apply to both continuous and transient operation of the sensors.

## 3 Corrections to the INTELLi-VENT-ASV Operator's Manual for HAMILTON-C1/T1

### 3.1 Vent Status and Quick Wean panels

In the **Vent Status** and **Quick Wean** panels, green bars are shown above and below the target zone when *all* the parameters are within the weaning zone.

When individual parameters are within their weaning zone, a timer is shown indicating the amount of time the parameter has been in the weaning zone.


### 3.2 Waveforms, view 7

When a CO2 sensor is connected, the following waveforms are shown in view 7:  
**Paw, Flow, PCO2**

When no CO2 sensor is connected the Volume waveform is shown in place of PCO2.

Note that the order of the waveforms on the display has changed (Figure 1).

Figure 1. Waveforms, view 7 in INTELLIVENT-ASV

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Hamilton Medical AG