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The screenshot shows the Scaffold DIA software interface. The main window displays a table of identified proteins. The table has columns for: #, Visible, Star, Protein Name, Molecular Weight, Protein Group Score, Identified Peptide Count, Exclusivity, Next CL Treatment, Taxonomy, Biological Process, Cellular Component, Molecular Function, Control, and Insulin. A 'Samples' button is highlighted in the left sidebar. The table lists various proteins such as HUMAN E3 SUMO-protein ligase, HUMAN Desmoglein OS=Homo sa..., HUMAN Alpha-actinin-4 OS=Hom..., HUMAN DNA topoisomerase 2-alp..., HUMAN Vinculin OS=Homo sapiens..., HUMAN Alpha-actinin-1 OS=Hom..., HUMAN Nucleolar RNA helicase 2..., HUMAN Ribosome-binding protei..., HUMAN Kinectin OS=Homo sapien..., HUMAN Splicing factor 3B subunit..., HUMAN Heat shock 70 kDa protein..., HUMAN Threonine-tRNA ligase, cycl..., HUMAN Eukaryotic translation initi..., HUMAN Trifunctional purine biosyn..., HUMAN T-complex protein 1 subun..., HUMAN T-complex protein 1 subun..., HUMAN Eukaryotic translation initi..., HUMAN Endoplasmic OS=Homo sa..., HUMAN Keratin, type II cytoskeletal..., HUMAN Phosphoglycerate kinase 1..., HUMAN Alpha-enolase OS=Homo..., and HUMAN Telomere-associated protei...

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
Create a chromatogram library and search against it : 2段階検索

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**Library Creation Parameters**

The following parameters were used in the conversion process and thus should match your instrumentation settings. If they vary dramatically from your settings please contact us.

Parameter	Setting
Charge Range	2 - 3
Maximum Missed Cleavages	1
m/z Range	396.4 - 1002.7
Default NCE	33
Default Charge	3

**Featured Libraries**

*Coronavirus reference*  
Download Coronavirus only DLIB (1.9 MB)  
Download Coronavirus plus Human (pan human) DLIB (259 MB)  
Download Coronavirus only FASTA - 13 entries  
Download Coronavirus plus Human FASTA - 20,350 entries

**Available Libraries**

*Arabidopsis thaliana*  
Download Arabidopsis thaliana DLIB (1.2 GB)  
Download FASTA file accessed 10/22/19 - 15,896 entries

*Caenorhabditis elegans*  
Download Caenorhabditis elegans DLIB (360 MB)  
Download FASTA file accessed 10/23/19 - 4,089 entries

*Danio rerio*  
Download Danio rerio DLIB (250 MB)  
Download FASTA file accessed 10/22/19 - 3,125 entries

*Drosophila melanogaster*  
Download Drosophila melanogaster DLIB (365 MB)  
Download FASTA file accessed 10/22/19 - 3,586 entries

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
## The Combined Workflow


For the highest quality identifications, Scaffold DIA can create a library. A pooled sample is searched with XCorDIA and the resulting ELIB is used to extract and quantify peptides from experimental samples.

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
### • Exclusive/total Intensity (Log<sub>2</sub> Fold Change)

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
### • Exclusive/total Quantified Peptide Count

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
Scaffold DIA - Demo2\_HeLa\_insulin\_6\_files\_30\_proteins.sdia

File Edit View Experiment Export Tools Help

Summarization: Treatment Protein FDR: 1.0% FDR Min # Peptides: 2

Filters: Show Hidden Name/Accession p-value filter GO Term

Display Type: Exclusive Intensity Normalized Log Intensities Color Options...

Color Legend (Displayed Value): 6.36E8, 1.22E8, 2.33E7

Proteins: 0.0% FDR (attained) 28 Targets 0 Decoys

Peptides: 0.0% FDR (attained) 880 Targets 0 Decoys

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Protein Sequence | Protein Level Charts | Chromatograms | Fragment Intensities | Fragmentation Table

sp|Q5T4S7|UBR4\_HUMAN E3 ubiquitin-protein ligase UBR4 OS=Homo sapiens GN=UBR4 PE=1 SV=1

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Q H N L L S P P F G W A S G S Q D S N S R R A T T P L ...

**H F S S D A V P H P R** F Y C V L S P E A S E D D L N R L D S V A C D V L F S K L 96

V K Y D E L Y A A L T A L L A A G S Q L D T V R R K E N K N V T A L E A C A ...

Y Y F L I L W R I L G I L P P S K **TYI N Q L S M N S P E M S E C D I L H T L R** 1040

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756 out of 5183 (14.58%) amino acids identified with 23 modifications

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Quantified	Peptide Sequence	Quantified Matches	Fixed Modif...	Variable Mo...	Start	Stop	Protein Accessi...	Proba...	Mass
<input type="checkbox"/>	AAPPPPPPPPLESSPR	0 of 6			606	622	sp Q5T4S7 UBR4...	86%	1,702.904
<input checked="" type="checkbox"/>	AEHASSLLELASTTK	0 of 6			1065	1079	sp Q5T4S7 UBR4...	100%	1,556.805
<input type="checkbox"/>	ALGTLGTTNEK	0 of 6			4803	4814	sp Q5T4S7 UBR4...	99%	1,234.623
<input type="checkbox"/>	APSYIEIFGR	0 of 6			2364	2373	sp Q5T4S7 UBR4...	98%	1,151.598
<input checked="" type="checkbox"/>	AGQALSELHVEK	2 of 6			1828	1850	sp Q5T4S7 UBR4...	100%	1,452.757
<input type="checkbox"/>									485.677
<input checked="" type="checkbox"/>									520.342
<input checked="" type="checkbox"/>									521.247

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②

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Sequence	Modifications	Charge	Sample	Quant. Intensity	# Qua...	RT Start (min)	RT Center (min)	RT Stop (min)	Precursor MZ	Attribute
3ALQYDTLISLMEHL...		3	Control_1	2.72E5	5	68.91	69.16	69.36	874.455	control
3ALQYDTLISLMEHL...		3	Control_2	—	0	68.55	—	69.32	874.455	control
ASVVTASSGSALQYDTLISLMEHL...		3	Control_3	4.455E5	5	68.65	68.89	69.13	874.455	control
ASVVTASSGSALQYDTLISLMEHL...		3	Control_2	3.8E5	5	68.97	69.21	69.41	874.455	insulin

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Protein Sequence | Protein Level Charts | Chromatograms | Fragment Intensities | Fragmentation Table

sp|Q5T4S7|UBR4\_HUMAN E3 ubiquitin-protein ligase UBR4 OS=Homo sapiens GN=UBR4 PE=1 SV=1

Q H N L L S P P F G W A S G S Q D S N S R R A T T P L Y H G F K E V E E N W S K 920



**H F S S D A V P H P R** F Y C V L S P E A S E D D L N R L D S V A C D V L F S K L 960

V K Y D E L Y A A L T A L L A A G S Q L D T V R R K E N K N V T A L E A C A L Q 1000

Y Y F L I L W R I L G I L P P S K **T Y I N Q L S M N S P E M S E C D I L H T L R** 1040

756 out of 5183 (14.59%) amino acids identified with 23 modifications

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Valid Sequence Modifications Charge Sample Quant. Intensity # Qua... RT Start (min) RT Center (min) RT Stop (min) Precu... Attrib...

QGFNWEVSGAGEAK 2 20170430\_HeLa\_DIA\_control\_2... 6.96185 5 41.04 41.37 41.63 789.892 Control

QGFNWEVSGAGEAK 2 20170430\_HeLa\_DIA\_insulin\_1... 7.94665 5 41.53 41.87 42.05 789.892 Insulin

QGFNWEVSGAGEAK 2 20170430\_HeLa\_DIA\_insulin\_2... 8.57665 5 41.55 41.88 42.14 789.892 Rapa...

Protein Sequence | Quantitative Charts | Chromatograms | Fragmentation Table

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Scaffold DIA Viewer - Demo2\_HeLa\_insulin\_6\_files\_30\_proteins.sdia

File Edit View Experiment Export Help

Summarization Condition Protein FDR 1.0% FDR Min # Peptides 2

Filters Show Hidden Name/Accession p-value filter GO Term

Organize Samples Proteins **Visualize** Analysis Publish

Proteins 0.0% FDR (attained) 29 Targets 0 Deays

Peptides 0.0% FDR (attained) 1002 Targets 0 Deays

**Volcano plot**

Principal Component Analysis (Log<sub>2</sub> Intensity)

Significance Marker 0-Fold Change Marker Multiselect Action Zoom Reference control

Log<sub>2</sub> p-Value vs Log<sub>2</sub> Fold Change (Precursor Intensity) (insulin/control)

Significant Sig. w/o correction Not significant

**Scatter plot**

Quantitative Scatterplot

contr vs insulin 45° -line Regression Line Multiselect Action Zoom


Log<sub>2</sub> Intensity (insulin) vs Log<sub>2</sub> Intensity (control)

Biological Process 3D Values Legend

Cellular process (28) Growth (1) Enzyme regulator activity (2) Biological adhesion (4) Developmental process (13) Establishment of localization (14) Biological regulation (21)

● Growth ● Enzyme regulator activity ● Biological adhesion ● Developmental process ● Biological regulation ● Cellular process

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



**Analysis**

Publish

Proteins  
0.0% FDR (attained)  
29 Targets  
0 Decoys

Peptides  
0.0% FDR (attained)  
1002 Targets  
0 Decoys

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Experiment Methods | SQL Report

Search

- Search Library: demo2elib
- Processing Directory: #work
- Protein Sequence Database: uniprot-swissprot-humanfasta
- Perform RT alignment between reference: false
- Fragmentation: CID
- Precursor Tolerance: 10.0 ppm
- Fragment Tolerance: 10.0 ppm
- Library Fragment Tolerance: 10.0 ppm
- Peptide FDR Threshold: 0.01
- Data Acquisition Type: Overlapping DIA
- Digestion Enzyme: Trypsin
- Peptide Length: [6 - 30]
- Peptide Charge: [2 - 3]
- Max Missed Cleavages: 1
- Modifications: Carbamidomethylation C 57.0214635 Non-termin...

Analysis

- Shared Evidence Clustering: Perfect
- Target Protein FDR: 0.01
- Minimum Number of Peptides: 2
- Grouping Applied in Version: 1.0.0
- Thresholding Applied in Version: 1.0.0
- Clustering Applied in Version: 1.0.0
- Quantify on Exclusive Peptides: true

Advanced

- Precursor Window Size: Deduced from file
- Minimum Number of Quant Ions: 3
- Maximum Number of Quant Ions: 5

Version

- Scaffold DIA: 1.0.0
- Encyclopedia: 0.6.12
- ProteoWizard: 3.0.11748
- Percolator: 3.0.1nightly-13-655e4c7-dirty

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**ANALYSIS OVERVIEW**

DIA data were analyzed using Scaffold DIA (1.0.0).

**RAW DATA PROCESSING**

Raw data files were converted to mzML format using ProteoWizard (3.0.11748). Deconvolution of overlapping windows was performed.

**SPECTRAL LIBRARY SEARCH**

Analytic samples were aligned based on retention times and individually searched against demo2.lib with a peptide mass tolerance of 10.0 ppm and a fragment mass tolerance of 10.0 ppm. Fixed modifications considered were: Carbamidomethylation C. The digestion enzyme was assumed to be Trypsin with a maximum of 1 missed cleavage site(s) allowed. Only peptides with charges in the range [2 - 3] and length in the range [6 - 30] were considered. Peptides identified in each sample were filtered by Percolator (3.0.1nightly-13-655e4c7-dirty) to achieve a maximum FDR of 0.01. Individual search results were combined and peptide identifications were assigned posterior error probabilities and filtered to an FDR threshold of 0.01 by Percolator (3.0.1nightly-13-655e4c7-dirty).

**QUANTIFICATION**

Peptide quantification was performed by Encyclopedia (0.6.12). For each peptide, the 5 highest quality fragment ions were selected for quantitation. Only peptides exclusive to each protein or cluster were used for quantitation.

**CRITERIA FOR PROTEIN IDENTIFICATION**

Proteins that contained similar peptides and could not be differentiated based on MS/MS analysis were grouped to satisfy the principles of parsimony. Proteins with a minimum of 2 identified peptides were thresholded to achieve a protein FDR threshold of 1.0%.

**GO ANNOTATION**

Proteins were annotated with GO terms from: UniProt, InterPro, GO\_Central, Reactome, GOC, HPA, Ensembl, IntAct, Pfam, UniProt, NTLU\_S6, LIPEDb, FlyBase, BHF-UCL, HNFIC, IKG1, SYSCILIA\_CONET, CACAO, AgBase, PINC, ARUK-UCL, CAFA, MTEBASE, Alzheimers\_University\_of\_Toronto, WormBase, GOB, SynGO-UCL, DFLAT, SGD, dictyBase and SynGO

**CITATIONS**


**ProteoWizard**  
A cross-platform toolkit for mass spectrometry and proteomics.  
Chambers, M.C., MacLean, B., Burke, R., Amodei, D., Ruderman, D.L., Neumann, S., Gatto, L., Fischer, B., Pratt, B., Egertson, J., Hoff, K., Kessner, D., Toman, N., Shulman, N., Frewen, B., Baker, T.A., Brusniak, M.-Y., Paulse, C., Oresny, D., Flashner, L., Kani, K., Moulden, C., Seymour, S.L., Nuwaysir, L.M., Lefebvre, B., Kuhlmann, F., Roark, J., Rainer, P., Detlev, S., Hemenway, T., Huhmer, A., Langridge, J., Connolly, B., Chadiack, T., Holly, K., Eckels, J., Deutsch, E.W., Moritz, R.L., Katz, J.E., Agus, D.B., MacCoss, M., Tabb, D.L. & Mallick, P. Nature Biotechnology 30, 918-920 (2012) <http://www.nature.com/nbt/journal/v30/n10/full/nbt.2377.html>

**Percolator**  
Semi-supervised learning for peptide identification from shotgun proteomics datasets  
Fujiwara, T., Koehn, C., Grotzer, T., Watanabe, M., Wilson, S., Bradford, A., and Mochly, J. MolCell

Copy Text to Clipboard | Export Publish Report | Export Supplementary Data

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ANOVA / t-test P ≥2 Background×Concentrations

Permutation Test NP ≥2 Background×Concentrations

Mann-Whitney U Test NP Exactly 2 Background×Concentrations

Kruskal-Wallis Test NP ≥2 Background×Concentrations

None

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## Multiple Test Correction

Control FWER with Hochberg's step-up and Holm's step-down ▾  
Control FDR with standard Benjamini-Hochberg procedure  
Control FWER with Hochberg's step-up and Holm's step-down  
No correction

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
- **Appendix A.** Structure of Scaffold DIA files (\*.sdia)
- **Appendix B.** Computation of FDR in Scaffold DIA
- **Appendix C.** Summarization: Rolling up Values
- **Appendix D.** Missing Values
- **Appendix E.** Shared Evidence Clustering Algorithm
- **Appendix F.** Heat map clustering
- **Appendix G.** Techniques to Control the Family-wise Error Rate
- **Appendix H.** Using Principal Component Analysis
- **Appendix I.** How PCA is Performed in Scaffold DIA
- **Appendix J.** Description of Mouse Right Click Context Menu Commands

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